

Hog 1000/500

Moving Light Consoles

TBJLT072.doc



Date: 9th May 2000

Hog 1000/500 PCB updates

Hog 1000/500 consoles now have updated CPU Cards to accommodate different MIDI and Keyboard connectors. The new metal connectors are now mounted on the CPU tray with wires and plugs connecting to sockets on the CPU card.

All Hog 1000 consoles from serial number B00565, and Hog 500 consoles from serial number C00145, have the updated CPU Card COLDCPU6 and metal DIN connectors.

If replacing the CPU card in consoles prior to the serial numbers listed above, COLDCPU6 may be used but will require the addition of four PCB mount DIN connectors ZPG4050.